

Nutrisafe[®] 2

Safety Enteral Feeding System*

Because their safety is no small matter, we've developed a **lockable non-luer system for safe enteral nutrition.**



We care... AS MUCH AS YOU DO.

Nutrisafe[®] 2

Safety Enteral Feeding System*



Non-Luer Design
A True Safety Connection

CM Markings on
Feeding Tubes

Ease of Use
During Placement

Secure Locking System
Prevents Disconnections
During Feeding

Vented Caps
Clean, Secure
Gravity Feeding

Optional Accessories Available
Suitable for Unique
Feeding Requirements

*The Nutrisafe 2 System is indicated for nasogastric/oralgastric enteral feeding, incorporating safety connectors which reduce the risk of misconnections between feeding tubes and intravenous connectors. The Nutrisafe 2 System is contraindicated for patients with congenital abnormalities of the gastrointestinal tract.

DEHP-Free Feeding Tubes

Length Units/Box Purple Fr 04 Gray Fr 05 Green Fr 06 Blue Fr 08 Black Fr 10 White Fr 12

PVC Feeding Tubes

50 cm	50	361.042	361.052	361.062	361.082	361.102	
75 cm	25	363.042	363.052	363.062	363.082	363.102	363.122
125 cm	25	362.042	362.052	362.062	362.082	362.102	362.122
40 cm	25	364.042	364.052	364.062	364.082		

Polyurethane Feeding Tubes

50 cm	50	1361.042	1361.052	1361.062	1361.082	1361.102	
75 cm	25	1363.042	1363.052	1363.062	1363.082	1363.102	1363.122
125 cm	25	1362.042	1362.052	1362.062	1362.082	1362.102	1362.122
40 cm	25	1364.042	1364.052	1364.062	1364.082		

Silicone Feeding Tubes

50 cm	10	2331.042		2331.062	2331.082	2331.102	
125 cm	10	2332.042		2332.062	2332.082	2332.102	
125 cm	10			2395.062	2395.082		



Syringes

Code	Size	Quantity
1015.002	.5 ml	100/box
1015.012	1 ml	100/box
1015.022	2.5 ml	100/box
1015.052	5 ml	100/box
1015.102	10 ml	100/box
1015.202	20 ml	50/box
1015.202M	20 ml	50/box
1015.352M	35 ml	25/box
1015.602	60 ml	25/box
1015.602M	60 ml	25/box



Extension Set (DEHP-free)

Code	Length	Diameter ID/OD	Qty
368.152	150 cm	1.5 / 2.5 mm	25/box
368.102	100 cm	1.5 / 2.5 mm	25/box
368.1521 [^]	150 cm	1.5 / 2.5 mm	25/box
368.1021 [^]	100 cm	1.5 / 2.5 mm	25/box
369.152	150 cm	3.0 / 4.2 mm	25/box

[^] includes clamps



Vygon USA 103A Park Drive, Montgomeryville, PA 18936

phone: 800-473-5414 fax: 215-672-6740 web: www.vygonusa.com

email: customerservice@vygonusa.com