


NEONATOLOGY
Vascular access



nutri**line**® &
nutri**line twinflo**®
2 Fr Polyurethane Neonatal PICCs



Value Life

2 Fr polyurethane neonatal PICCs

nutriline®, 2 Fr peripherally inserted central venous catheter with integrated extension line and fixation wing, allows for long term venous access for premature babies, neonates and infants.

nutriline twinflo® double-lumen allows for infusion of incompatible drugs or infusion of medication and parenteral nutrition solutions at the same time through a single access point.

- Thermosensitive polyurethane material softens at body temperature.
- Radiopaque catheter provides clear visualization for terminal tip placement.
- Catheter markings every centimeter.
- Available with choice of introducers including splitting needle and microflash® peelable cannula.



Nutriline



Product Code	Size	Length	Priming volume	Flow rate	Introducer	Stylet	Units/Case
1252.30G	2 Fr	30 cm	0.12 ml	5.0 ml / min	Splitting needle 20 G	No	10
1252.31G	2 Fr	30 cm	0.12 ml	5.0 ml / min	No introducer	No	10
1252.35G	2 Fr	30 cm	0.12 ml	5.0 ml / min	Microflash 20 G	No	10
1252.030G	2 Fr	30 cm	0.12 ml	5.0 ml / min	Microflash 20 G	Yes	10

Nutriline Twinflo

Product Code	Catheter			Introducer		Stylet	Units/Case
	Length	Priming volume	Flow rate (pressure lbar)	Type	Size		
1252.232G	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	No introducer	-	No	10
1252.232M	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	Microflash	20 G	No	10
1252.230G	30 cm	2 x 0.2 ml	2 x 1.45 ml / min	Microflash	20 G	Yes	10

Kits & accessories

Product Code	Description	Units/Case
7371.19	20 G Microflash introducer	10
7371.195	20 G Microflash safety introducer	10
1262.30N	20 G Splitting needle	20
AMS-8198-I	Neonatal PICC insertion tray	20
AMS-8198CS	Neonatal PICC insertion tray with CHG	20
AMS-8199	Neonatal PICC insertion tray with V-cutter	20
5804.08	Grip-Lok® catheter securement device	20
AMS-850TC	V-cutter catheter trimming tool	20



For further information, please contact: customerservice@vygonusa.com.

The specifications shown in this brochure are for information only and are not, under any circumstances, of a contractual nature.