


NEONATOLOGY

Vascular access



premicath[®]
1Fr Polyurethane
Neonatal PICC





Value Life

1 Fr polyurethane neonatal PICC


premicath®, 1 Fr peripherally inserted central venous catheter with integrated extension line and fixation wing, allows for long term venous access ideal for premature and low-weight babies (less than 1kg) for delivery of medication and parenteral nutrition.

- Thermosensitive polyurethane material softens at body temperature.
- Radiopaque catheter provides clear visualization for terminal tip placement.
- Catheter markings every centimeter.
- Available with choice of introducers including splitting needle, microflash® peelable cannula and IV cannula.



Premicath

Product Code	Catheter			Introducer		Stylet	Units/Case
	Length	Priming volume	Flow rate (pressure lbar)	Type	Size		
1261.153G	15 cm	0.08 ml	0.9 ml / min	IV cannula	24 G	Yes	10
1261.203G	20 cm	0.09 ml	0.7 ml / min	Splitting needle	24 G	Yes	10
1261.207G	20 cm	0.09 ml	0.7 ml / min	IV cannula	24 G	Yes	10
1261.20G	20 cm	0.09 ml	0.7 ml / min	Splitting needle	24 G	No	10
1261.208G	20 cm	0.09 ml	0.7 ml / min	Microflash	20 G	Yes	10
1261.306G	30 cm	0.11 ml	0.6 ml / min	No introducer	-	Yes	10



Kits & accessories

Product Code	Description	Units/Case
7371.19	20 G Microflash introducer	10
7371.195	20 G Microflash safety introducer	10
1261.20N	24 G Splitting needle	20
CMS-207	24 G IV cannula	10
AMS-8198-I	Neonatal PICC insertion tray	20
AMS-8198CS	Neonatal PICC insertion tray with CHG	20
AMS-8199	Neonatal PICC insertion tray with V-cutter	20
5804.08	Grip-Lok® catheter securement device	20
AMS-850TC	V-cutter catheter trimming tool	20

For further information, please contact: customerservice@vygonusa.com.

The specifications shown in this brochure are for information only and are not, under any circumstances, of a contractual nature.