


NEONATOLOGY
Vascular access

Umbilical Catheters Immediate Central Access





Value Life

Polyurethane umbilical catheters

For babies requiring catheterization, polyurethane umbilical catheters allow for immediate central access with easy and comfortable insertion.

- Thermosensitive polyurethane material softens at body temperature.
- Radiopaque catheter provides clear visualization for terminal tip placement.
- Double-lumen with staggered exit points ideal for infusing incompatible drugs.



Single-lumen catheters

Product Code	Size	Length	Int. ø	Priming Volume	Flow rate	Units/Case
1270.02	2.5 Fr	30 cm	0.5 mm	0.21 ml	2 ml / min	8
1270.03	3.5 Fr	40 cm	0.8 mm	0.34 ml	11 ml / min	8
1270.04	4 Fr	40 cm	0.8 mm	0.36 ml	18 ml / min	20
1270.05	5 Fr	40 cm	1.0 mm	0.46 ml	24 ml / min	20
1270.08	8 Fr	40 cm	1.5 mm	0.84 ml	79 ml / min	20

Double-lumen catheters

Product Code	Size	Length	Lumen	Priming Volume		Flow rate		Units/Case
				Distal	Proximal	Distal	Proximal	
1272.14	4 Fr	20 cm	2 x 20G	0.26 ml	0.26 ml	15 ml / min	15 ml / min	10
1274.14	4 Fr	40 cm	2 x 20G	0.28 ml	0.28 ml	6 ml / min	6 ml / min	10
1274.17	5 Fr	40 cm	2 x 19G	0.30 ml	0.30 ml	9 ml / min	9 ml / min	10

Kits & accessories

Product Code	Description	Units/Case
AMS-9294	Umbilical catheter insertion tray	10
AMS-9294CS	Umbilical catheter insertion tray with CHG	10

For further information, please contact: customerservice@vygonusa.com.

The specifications shown in this brochure are for information only and are not, under any circumstances, of a contractual nature.